

**Certificate of Analysis
GM501 Flush**

Batchnumber	FLUSH-2206
Production date	17 June 2022
Expiry date	31 December 2022
Total volume (mL)	71600

	Results	Specifications
pH	7,39	7,20-7,40
Osmolality	279	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date Dr. Julia Heinzmann

Released by 07 July 2022

Signature

p.p. Julia Heinzmann