

Certificate of Analysis GM501 Flush

Batchnumber	FLUSH-2209
Production date	20 October 2022
Expiry date	30 April 2023
Total volume (mL)	88900

	Results	Specifications
pH	7,36	7,20-7,40
Osmolality	279	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date	04 November 2022
Released by	Nicole Geiger
Signature	