

**Certificate of Analysis
GM501 Flush**


Batchnumber	FLUSH-2211
Production date	06 December 2022
Expiry date	30 June 2023
Total volume (mL)	77700

	Results	Specifications
pH	7,36	7,20-7,40
Osmolality	280	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date 21 December 2022

Released by Nicole Geiger

Signature



pp. N. Geiger