

**Certificate of Analysis
GM501 Flush**

Batchnumber	FLUSH-2212
Production date	14 December 2022
Expiry date	30 June 2023
Total volume (mL)	51200

	Results	Specifications
pH	7,32	7,20-7,40
Osmolality	285	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date 02 January 2023

Released by Nicole Geiger

Signature

