

**Certificate of Analysis
GM501 Flush**

Batchnumber	FLUSH-2307
Production date	04 July 2023
Expiry date	31 January 2024
Total volume (mL)	67600

	Results	Specifications
pH	7,37	7,20-7,40
Osmolality	281	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date	19 July 2023
--------------	--------------

Released by	Dr. Julia Heinzmann
-------------	---------------------

Signature

p.p. Julia Heinzmann