

**Certificate of Analysis
GM501 Flush**

Batchnumber	FLUSH-2308
Production date	10 August 2023
Expiry date	28 February 2024
Total volume (mL)	87500

	Results	Specifications
pH	7,37	7,20-7,40
Osmolality	283	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	90,48	≥ 80 %

Release date	28 August 2023
Released by	Nicole Geiger
Signature	

pp. N. Geiger