

**Certificate of Analysis
GM501 Flush**

Batchnumber FLUSH-2309
 Production date 27 September 2023
 Expiry date 31 March 2024
 Total volume (ml) 86500

Parameter	Results	Specifications
pH	7,35	7,20-7,40
Osmolality	282	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse embryo test (1-cell)	90,48	≥ 80 %

Release date 12 October 2023
 Released by Nicole Geiger
 Signature

pp. N. Geiger

REV 07_00 / 2023-09-15