

**Certificate of Analysis
GM501 Flush**

Batchnumber FLUSH-2310
 Production date 14 November 2023
 Expiry date 31 May 2024
 Total volume (ml) 108600

Parameter	Results	Specifications
pH	7,32	7,20-7,40
Osmolality	280	270-290 mOsm/kg
Bacterial screen	No growth	No growth
Fungal screen	No growth	No growth
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse embryo test (1-cell)	90,48	≥ 80 %

Release date 29 November 2023
 Released by Nicole Geiger
 Signature

pp. N. Geiger

REV 07_00 / 2023-09-15