

## Certificate of Analysis GM501 Flush

Batchnumber	FLUSH-2404
Production date	14 March 2024
Expiry date	30 September 2024
Total volume (ml)	13200

Parameter	Results	Specifications
pH	7,24	7,20-7,40
Osmolality	280	270-290 mOsm/kg
Bacterial screen	NO GROWTH	No growth
Fungal screen	NO GROWTH	No growth
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse embryo test <small>(1-cell)</small>	100	≥ 80 %

Release date	02 April 2024
Released by	Nicole Geiger
Signature	

REV 07\_00 / 2023-09-15