

**Certificate of Analysis
GM501 Flush**

Batchnumber FLUSH-2405
 Production date 08 April 2024
 Expiry date 31 October 2024
 Total volume (ml) 75300

Parameter	Results	Specifications
pH	7,31	7,20-7,40
Osmolality	282	270-290 mOsm/kg
Bacterial screen	No growth	No growth
Fungal screen	No growth	No growth
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse embryo test (1-cell)	90,48	≥ 80 %

Release date 30 April 2024
 Released by Joana Sturm
 Signature



REV 07_00 / 2023-09-15