


**Certificate of Analysis
GM501 Flush**

Batchnumber FLUSH-2406
 Production date 02 May 2024
 Expiry date 30 November 2024
 Total volume (ml) 81600

Parameter	Results	Specifications
pH	7,33	7,20-7,40
Osmolality	278	270-290 mOsm/kg
Bacterial screen	NO GROWTH	No growth
Fungal screen	NO GROWTH	No growth
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse embryo test (1-cell)	86	≥ 80 %

Release date 17 May 2024
 Released by Nicole Geiger
 Signature



REV 07_00 / 2023-09-15