


**Certificate of Analysis  
GM501 Flush**

Batchnumber FLUSH-2407  
 Production date 06 June 2024  
 Expiry date 31 December 2024  
 Total volume (ml) 93900

Parameter	Results	Specifications
pH	7,31	7,20-7,40
Osmolality	277	270-290 mOsm/kg
Bacterial screen	No growth	No growth
Fungal screen	No growth	No growth
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse embryo test (1-cell)	100	≥ 80 %

Release date 25 June 2024  
 Released by Nicole Geiger  
 Signature



REV 07\_00 / 2023-09-15