

**Certificate of Analysis
GM501 Flush**

Batchnumber FLUSH-2408
 Production date 02 July 2024
 Expiry date 31 January 2025
 Total volume (ml) 54500

Parameter	Results	Specifications
pH	7,3	7,20-7,40
Osmolality	279	270-290 mOsm/kg
Bacterial screen	No growth	No growth
Fungal screen	No growth	No growth
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse embryo test (1-cell)	100	≥ 80 %

Release date 22 July 2024
 Released by Nicole Geiger
 Signature

pp. N. Geiger

REV 07_00 / 2023-09-15