

**Certificate of Analysis  
GM501 Flush**

Batchnumber FLUSH-2409  
 Production date 02 August 2024  
 Expiry date 28 February 2025  
 Total volume (ml) 75400

Parameter	Results	Specifications
pH	7,3	7,20-7,40
Osmolality	281	270-290 mOsm/kg
Sterility test (performed according to the current PH. Eur. 2.6.1 / USP <71>)	No growth	No growth
Concentration heparin	2,68	2-3,75 IU/mL
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse embryo test (1-cell)	90,48	≥ 80 %

Release date 19 August 2024

Released by Joana Sturm

Signature



REV 08\_00 / 2024-08-16