

**Certificate of Analysis
GM501 Flush**

Batchnumber FLUSH-2412
 Production date 02 October 2024
 Expiry date 30 April 2025
 Total volume (ml) 84400

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH	7,26	7,20-7,40
Osmolality	280	270-290 mOsm/kg
Sterility test (performed according to the current PH. Eur. 2.6.1 / USP <71>)	No growth	No growth
Concentration heparin	2,61	2-3,75 IU/mL
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse embryo test (1-cell)	95,24	≥ 80 %

Release date 21 October 2024

Released by Joana Sturm

Signature



REV 08_00 / 2024-08-16