

**Certificate of Analysis
GM501 Flush**

Batchnumber FLUSH-2414
 Production date 13 December 2024
 Expiry date 30 June 2025
 Total volume (ml) 132.400

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH	7,33	7,20-7,40
Osmolality	278	270-290 mOsm/kg
Sterility test (performed according to the current PH. Eur. 2.6.1 / USP <71>)	No growth	No growth
Concentration heparin	2,28	2-3,75 IU/mL
Endotoxin level	< 0,05	< 0,25 EU/mL
MEA (1-cell to expanded blastocyst within 96 hours)	90,48	≥ 80 %

Release date 02 January 2025

Released by Nicole Geiger

Signature



REV 08_00 / 2024-08-16