


Certificate of Analysis GM501 Flush

Batchnumber	FLUSH-2505
Production date	02 June 2025
Expiry date	31 December 2025
Total volume (ml)	93000
Product was sterilized by aseptic processing techniques	

Parameter	Results	Specifications
pH	7,31	7,20-7,40
Osmolality	283	270-290 mOsm/kg
Sterility test (performed according to the current PH. Eur. 2.6.1 / USP <71>)	No growth	No growth
Concentration heparin	2,22	2-3,75 IU/mL
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse embryo test (1-cell)	100	≥ 80 %

Release date	23 June 2025
Released by	Nicole Geiger
Signature	

REV 08_00 / 2024-08-16