

Certificate of Analysis GM501 Flush

Batchnumber FLUSH-2605

Production date 13 April 2026

Expiry date 31 October 2026

Total volume (ml) 90200

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH	7,30	7,20-7,40
Osmolality	281	270-290 mOsm/kg
Sterility test (performed according to the current PH. Eur. 2.6.1 / USP <71>)	No growth	No growth
Concentration heparin	2,48	2-3,75 IU/mL
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse embryo test (1-cell)	95	≥ 80 %

Release date 29 April 2026

Released by Dr. René Richarz

Signature



REV 08_00 / 2024-08-16