

Certificate of Analysis
GM501 Gradient 45%

Batchnumber	GRAD45-2353
Production date	09 November 2023
Expiry date	31 May 2025
Total volume (ml)	16450

Parameter	Results	Specifications
pH	7,37	7,20-7,60
Osmolality	329	310-340 mOsm/kg
Endotoxin level	< 0,1	< 0,50 EU/mL
Sterility test (according to the current Ph. Eur. 2.6.1 / <71>)	No growth	No growth
Sperm survival Test (4 hours exposure)	98	≥ 80 %
Sperm survival Test (24 hours exposure)	86	≥ 75 %

Release date	27 November 2023
Released by	Nicole Geiger
Signature	

pp. N. Geiger

REV 08_00 / 2023-09-15