


**Certificate of Analysis
GM501 Gradient 45%**

Batchnumber GRAD45-2409
 Production date 04 April 2024
 Expiry date 31 October 2025
 Total volume (ml) 9150

Parameter	Results	Specifications
pH	7,39	7,20-7,60
Osmolality	323	310-340 mOsm/kg
Endotoxin level	< 0,15	< 0,50 EU/mL
Sterility test (according to the current Ph. Eur. 2.6.1 / <71>)	NO GROWTH	No growth
Sperm survival Test (4 hours exposure)	96	≥ 80 %
Sperm survival Test (24 hours exposure)	88	≥ 75 %

Release date 25 April 2024
 Released by Nicole Geiger

Signature



REV 08_00 / 2023-09-15