

Certificate of Analysis

GM501

Gradient 90%


Batchnumber GRAD90-2416
 Production date 02 July 2024
 Expiry date 31 January 2026
 Total volume (ml) 13250

Parameter	Results	Specifications
pH	7,38	7,20-7,60
Osmolality	336	320-350 mOsm/kg
Density	1,1103	1,1050-1,1150 g/ml
Endotoxin level	< 0,15	< 0,50 EU/ml
Sterility test (according to the current Ph. Eur. 2.6.1 / USP <71>)	No growth	No growth
Viscosity in cP (25°C)	1,54	< 1,65 cP
Sperm survival Test (4 hours exposure)	100	≥ 80 %
Sperm survival Test (4 hours exposure)	86	≥ 75 %

Release date 17 July 2024

Released by Nicole Geiger

Signature



REV 08_00 / 2023-09-15