

**Certificate of Analysis  
GM501 Gradient 100%**

Batchnumber	GSS-2114
Production date	11 June 2021
Expiry date	31 December 2022
Total volume (mL)	80000

	<b>Results</b>	<b>Specifications</b>
pH	7,35	7,20-7,60
Osmolality	321	300-330 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,1	< 0,50 EU/mL
Density	1,1212	1,1150-1,1250 g/mL
Viscosity in cP (26°C)	1,73	< 1,75 cP
Sperm survival Test (4 hours exposure)	98	≥ 80 %
Sperm survival Test (24 hours exposure)	95	≥ 75 %

Release date 28 June 2021

Released by Vanessa Bierlich

Signature