

Certificate of Analysis
GM501 Gradient 100%

Batchnumber GSS-2340
Production date 08 November 2023
Expiry date 31 May 2025
Total volume (ml) 209750

Parameter	Results	Specifications
pH	7,34	7,20-7,60
Osmolality	325	300-330 mOsm/kg
Endotoxin level	< 0,1	< 0,50 EU/ml
Density	1,1197	1,1150-1,1250 g/ml
Sterility test (according to the current Ph.Eur. 2.6.1 / USP <71>)	No growth	No growth
Viscosity in cP (25°C)	1,64	< 1,75 cP
Sperm survival Test (4 hours exposure)	97	≥ 80 %
Sperm survival Test (24 hours exposure)	88	≥ 75 %

Release date 23 November 2023

Released by Nicole Geiger

Signature

pp. N. Geiger

REV 07_00 / 2023-09-15