

**Certificate of Analysis
GM501 Gradient 100%**

Batchnumber GSS-2619
Production date 26 June 2026
Expiry date 31 December 2027
Total volume (ml) 78110

Product is sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH	7,29	7,20-7,60
Osmolality	314	300-330 mOsm/kg
Endotoxin level	< 0,15	< 0,50 EU/ml
Density	1,1189	1,1150-1,1250 g/ml
Sterility test (according to the current Ph.Eur. 2.6.1 / USP <71>)	No growth	No growth
HSSA (% motility compared with control after 4h exposure to test medium)	95	≥ 80 %
HSSA (% motility compared with control after 24h exposure to test medium)	90	≥ 75 %

Release date 14 July 2026

Released by Nicole Geiger

Signature

pp. N. Geiger

REV 08_00 / 2025-01-10