

**Certificate of Analysis
GM501 HSA**

Batchnumber	HSA-2201
Production date	08 March 2022
Expiry date	31 March 2023
Total volume (mL)	10210
Batchnumber HSA	K151A6872

	Results	Specifications
pH	7,3	7,20-7,40
Osmolality	279	265-285 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,8	< 1,0 EU/mL
Color	AMBER	AMBER
Clarity	Clear, no particles	Clear, no particles
Mouse Embryo Assay (1-cell)	95	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	30 March 2022
Released by	Dr. Julia Heinzmann
Signature	

p.p. Julia Heinzmann