

**Certificate of Analysis  
GM501 HSA**

Batchnumber	HSA-2301
Production date	24 January 2023
Expiry date	31 January 2024
Total volume (mL)	7000
Batchnumber HSA	K235A687H

	<b>Results</b>	<b>Specifications</b>
pH	7,35	7,20-7,40
Osmolality	279	265-285 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,8	< 1,0 EU/mL
Color	AMBER	AMBER
Clarity	PASSED	PASSED
Mouse Embryo Assay (1-cell)	100	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date	08 February 2023
Released by	Nicole Geiger
Signature	

REV 05\_00 /2019-01-15