

## Certificate of Analysis GM501 HSA

Batchnumber	HSA-2403
Production date	15 April 2024
Expiry date	30 April 2025
Total volume (ml)	10500
Batchnumber HSA	K331A687A

Parameter	Results	Specifications
pH	7,26	7,20-7,40
Osmolality	276	265-285 mOsm/kg
Endotoxin	< 0,2	< 1,0 EU/mL
Color	AMBER	AMBER
Sterility (fungal and bacterial sterility test)	PASSED	PASSED
Clarity	Clear, no particles	Clear, no particles
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	100	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date	30 May 2024
Released by	Nicole Geiger
Signature	

REV 06\_00 / 2023-09-15