

Certificate of Analysis
GM501 HSA

Batchnumber	HSA-2501
Production date	03 March 2025
Expiry date	31 March 2026
Total volume (ml)	10000
Batchnumber HSA	M424A6871

Parameter	Results	Specifications
pH	7,35	7,20-7,40
Osmolality	284	265-285 mOsm/kg
Endotoxin	< 0,2	< 1,0 EU/mL
Color	AMBER	AMBER
Sterility (fungal and bacterial sterility test)	PASSED	PASSED
Clarity	Clear, no particles	Clear, no particles
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	19 March 2025
Released by	Nicole Geiger
Signature	

pp. N. Geiger

REV 06_00 / 2023-09-15