

Certificate of Analysis GM501 HSA

Batchnumber	HSA-2603
Production date	04 May 2026
Expiry date	31 May 2027
Total volume (ml)	8250
Batchnumber HSA	KP28A6871

Parameter	Results	Specifications
pH	7,28	7,20-7,40
Osmolality	284	265-285 mOsm/kg
Endotoxin level	< 0,3	< 1,0 EU/mL
Concentration HSA	98	80-120 g/l
Color	Amber	AMBER
Sterility (fungal and bacterial sterility test)	No growth	NO GROWTH
Clarity	Clear, no particles	Clear, no particles
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	20 May 2026
Released by	Nicole Geiger
Signature	

REV 07_00 / 2025-10-01