

**Certificate of Analysis
GM501 HSA**

Batchnumber HSA-2402
 Production date 22 February 2024
 Expiry date 28 February 2025
 Total volume (ml) 9000
 Batchnumber HSA M306B687M

Parameter	Results	Specifications
pH	7,26	7,20-7,40
Osmolality	277	265-285 mOsm/kg
Endotoxin	< 0,85	< 1,0 EU/mL
Color	AMBER	AMBER
Sterility (fungal and bacterial sterility test)	PASSED	PASSED
Clarity	Clear, no particles	Clear, no particles
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 25 April 2024
 Released by Nicole Geiger
 Signature

