

**Certificate of Analysis
GM501 Hyaluronidase**

Batchnumber HYA-1928
 Production date 05 December 2019
 Expiry date 31 December 2020
 Total volume (ml) 3600
 Batchnumber HSA B3ALC00833 / B3ALD00533

	Results	Specifications
pH	7,43	7,30-7,60
Osmolality	280	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,112	< 1,00 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 02 January 2020

Released by Vanessa Bierlich

Signature