

**Certificate of Analysis  
GM501 Hyaluronidase**

Batchnumber HYA-2020  
 Production date 25 November 2020  
 Expiry date 30 November 2021  
 Total volume (ml) 3600  
 Batchnumber HSA B3ALE00033

	<b>Results</b>	<b>Specifications</b>
pH	7,45	7,30-7,60
Osmolality	279	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,184	< 1,00 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 14 December 2020

Released by Vanessa Bierlich

Signature