

**Certificate of Analysis
GM501 Hyaluronidase**

Batchnumber	HYA-2114
Production date	21 June 2021
Expiry date	30 June 2022
Total volume (ml)	3600
Batchnumber HSA	K041A6871/K052B6871

	Results	Specifications
pH	7,44	7,30-7,60
Osmolality	282	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,149	< 1,00 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 06 July 2021

Released by Vanessa Bierlich

Signature