

**Certificate of Analysis
GM501 Hyaluronidase**

Batchnumber	HYA-2218
Production date	12 August 2022
Expiry date	31 August 2023
Total volume (ml)	3700
Batchnumber HSA	B3ALF01143
Batchnumber HYA	S1026H40

	Results	Specifications
pH	7,49	7,30-7,60
Osmolality	281	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,179	< 1,00 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 30 August 2022

Released by Dr. Julia Heinzmann

Signature

p.p. Julia Heinzmann

REV 08_00 /2022-05-09