

**Certificate of Analysis  
GM501 Hyaluronidase**

Batchnumber	HYA-2228
Production date	20 December 2022
Expiry date	31 December 2023
Total volume (ml)	3800
Batchnumber HSA	K222B6875
Batchnumber HYA	S1026H40

	<b>Results</b>	<b>Specifications</b>
pH	7,52	7,30-7,60
Osmolality	285	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,2	< 1,00 EU/mL
Mouse Embryo Assay (1-cell)	95	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 05 January 2023

Released by Nicole Geiger

Signature

