

**Certificate of Analysis**  
**GM501 Hyaluronidase**

Batchnumber	HYA-2326
Production date	26 September 2023
Expiry date	30 September 2024
Total volume (ml)	3600
Batchnumber HSA	B03H003853
Batchnumber HYA	S1047H40

Parameter	Results	Specifications
pH	7,47	7,30-7,60
Osmolality	279	270-290 mOsm/kg
Endotoxin	< 0,1	<0,5 EU/ml
Concentration HSA	3,9	3,2-4,8 g/l
Enzyme activity	75-110	40-160 IU/ml
Sterility test (according to the current Ph.Eur.2.6.1 /USP<71>)	No growth	No growth
Mouse Embryo Assay (1-cell)	80,95	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	23 October 2023
Released by	Nicole Geiger
Signature	

REV 09\_00 /2023-09-15