


**Certificate of Analysis
GM501 Hyaluronidase**

Batchnumber	HYA-2624
Production date	29 June 2026
Expiry date	30 June 2027
Total volume (ml)	6000
	B23K049493
Batchnumber HSA	B23K064693
Batchnumber HYA	S2015H40-300320
Product was sterilized by aseptic processing techniques	

Parameter	Results	Specifications
pH	7,49	7,30-7,60
Osmolality	276	270-290 mOsm/kg
Endotoxin	< 0,13	<0,5 EU/ml
Concentration HSA	3,9	3,2-4,8 g/l
Enzyme activity	56-83	40-160 IU/ml
Sterility test (according to the current Ph.Eur.2.6.1 /USP<71>)	No growth	No growth
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	15 July 2026
Released by	 Nicole Geiger
Signature	

REV 09_00 /2023-09-15