

## Certificate of Analysis GM501 Mineral Oil

Batchnumber	OIL-1704
Production date	30 January 2017
Expiry date	31 July 2018
Total volume (mL)	35000

	Results	Specifications
Viscosity	20,3	< 30 cP at 30°C
Appearance	PASSED	PASSED
Mouse Embryo Assay <small>(1-cell to expanded blastocyst in 96 hours)</small>	95,24	≥ 80 %
Endotoxin	< 0,1	< 0,25 EU/mL
Sterility test <small>Sterilised by sterile filtration over 0,2µm filter</small>	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Specific gravity	0,8446	< 0,86 g/mL
POV	< 0,1	< 0,1 mEQ/Kg

Release date 14 February 2017

Released by Dr. Julia Heinzmann

Signature

*p.p. Julia Heinzmann*