

**Certificate of Analysis
GM501 Mineral Oil**

Batchnumber	OIL-1709
Production date	15 March 2017
Expiry date	30 September 2018
Total volume (ml)	50000

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	20,6	< 30 cP at 30°C
Specific gravity	0,8458	< 0,86 g/ml
Endotoxin	< 0,1	< 0,25 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	95,24	≥ 80 %

Release date 31 March 2017

Released by Dr. Julia Heinzmann

Signature

p.p. Julia Heinzmann