

**Certificate of Analysis  
GM501 Mineral Oil**

Batchnumber	OIL-1724
Production date	07 August 2017
Expiry date	28 February 2019
Total volume (ml)	33200

	<b>Results</b>	<b>Specifications</b>
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	22,1	< 30 cP at 30°C
Specific gravity	0,8476	< 0,86 g/ml
Endotoxin	< 0,1	< 0,25 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %

Release date 23 August 2017

Released by Dr. Julia Heinzmann

Signature

*p.p. Julia Heinzmann*