

**Certificate of Analysis
GM501 Mineral Oil**

Batchnumber OIL-1726
 Production date 29 August 2017
 Expiry date 28 February 2019
 Total volume (ml) 35000

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	22,3	< 30 cP at 30°C
Specific gravity	0,8479	< 0,86 g/ml
Endotoxin	< 0,1	< 0,25 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %

Release date 15 September 2017

Released by Vanessa Bierlich

Signature