

Certificate of Analysis GM501 Mineral Oil

Batchnumber	OIL-1805
Production date	05 February 2018
Expiry date	31 August 2019
Total volume (ml)	15000

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	22,5	< 30 cP at 30°C
Specific gravity	0,8494	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	90,48	≥ 80 %

Release date	21 February 2018
Released by	Dr. Julia Heinzmann
Signature	

p.p. Julia Heinzmann