

**Certificate of Analysis  
GM501 Mineral Oil**

Batchnumber	OIL-1826
Production date	06 September 2018
Expiry date	31 March 2020
Total volume (ml)	199500

	<b>Results</b>	<b>Specifications</b>
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	21,9	< 30 cP at 30°C
Specific gravity	0,8435	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	90,48	≥ 80 %

Release date	21 September 2018
Released by	Dr. Julia Heinzmann
Signature	

*p.p. Julia Heinzmann*