

**Certificate of Analysis  
GM501 Mineral Oil**

Batchnumber	OIL-1831
Production date	25 October 2018
Expiry date	30 April 2020
Total volume (ml)	190000

	<b>Results</b>	<b>Specifications</b>
Appearance	PASSED	PASSED
Sterility test <small>Sterilised by sterile filtration over 0,2µm filter</small>	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	20,8	< 30 cP at 30°C
Specific gravity	0,8428	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay <small>(1-cell to expanded blastocyst after 120 hours)</small>	80,95	≥ 80 %

Release date 12 November 2018

Released by Dr. Julia Heinzmann

Signature

*p.p. Julia Heinzmann*