

**Certificate of Analysis  
GM501 Mineral Oil**

Batchnumber	OIL-1901
Production date	09 January 2019
Expiry date	31 July 2020
Total volume (mL)	179000

	<b>Results</b>	<b>Specifications</b>
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	24,5	< 30 cP at 30°C
Specific gravity	0,852	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %

Release date 24 January 2019

Released by Vanessa Bierlich

Signature