

**Certificate of Analysis
GM501 Mineral Oil**

Batchnumber	OIL-1921
Production date	18 June 2019
Expiry date	31 December 2020
Total volume (mL)	182000

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	20,9	< 30 cP at 30°C
Specific gravity	0,8487	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	85,71	≥ 80 %

Release date	03 July 2019
Released by	Vanessa Bierlich

Signature

