

**Certificate of Analysis  
GM501 Mineral Oil**

Batchnumber OIL-1939  
 Production date 20 November 2019  
 Expiry date 31 May 2021  
 Total volume (mL)

	<b>Results</b>	<b>Specifications</b>
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	22	< 30 cP at 30°C
Specific gravity	0,8429	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	90,48	≥ 80 %

Release date 05 December 2019  
 Released by Vanessa Bierlich

Signature