

**Certificate of Analysis  
GM501 Mineral Oil**

Batchnumber	OIL-2006
Production date	11 February 2020
Expiry date	31 August 2021
Total volume (mL)	168500

	<b>Results</b>	<b>Specifications</b>
Appearance	PASSED	PASSED
Sterility test <small>Sterilised by sterile filtration over 0,2µm filter</small>	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	20,8	< 30 cP at 30°C
Specific gravity	0,8429	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay <small>(1-cell to expanded blastocyst after 120 hours)</small>	100	≥ 80 %

Release date	02 March 2020
Released by	Vanessa Bierlich

Signature