

Certificate of Analysis GM501 Mineral Oil

Batchnumber	OIL-2016
Production date	29 June 2020
Expiry date	31 December 2021
Total volume (mL)	169100

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test <small>Sterilised by sterile filtration over 0,2µm filter</small>	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	22	< 30 cP at 30°C
Specific gravity	0,8429	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay <small>(1-cell to expanded blastocyst after 120 hours)</small>	95,83	≥ 80 %

Release date 14 July 2020

Released by Dr. Julia Heinzmann

Signature

p.p. Julia Heinzmann