

**Certificate of Analysis
GM501 Mineral Oil**

Batchnumber	OIL-2037
Production date	08 December 2020
Expiry date	30 June 2022
Total volume (mL)	185000

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	20,5	< 30 cP at 30°C
Specific gravity	0,8408	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	85,95	≥ 80 %

Release date	23 December 2020
Released by	Vanessa Bierlich

Signature